

# Laboratories Swab Identification Guide Health New Zealand



#### **Bacterial Swabs**

For stock orders email labinfo@cdhb.health.nz. For other enquires phone Microbiology ext. 80350



Bacterial Swabs (purple) + Transport Medium for Routine Culture (including cervical swabs for N. gonorrhoeae)

CERVICAL SWABS / BACTERIAL SWABS (as above) for Mycoplasma hominis / Ureaplasma urealyticum. Place into specialised Mycoplasma/Ureaplasma medium.



RILE CONTAI

#### Chlamydia trachomatis, N. gonorrhoea, T. vaginalis, M. genitalium

Alinity-m multi-Collect specimen collection kit. Use swab for all genital sites, rectum and eyes.

Use transfer pipette for aliquoting urine up to window on multi-Collect Tube.

Refer "specimen collection guides" how to sample for vaginal, urine, rectal specimens.

**CAUTION:** Use of lubricant gels can cause PCR inhibition. If gels cannot be avoided, then please use a non-interfering

gel such as "Shield Gel".



### Mycoplasma.genitalium self-collection kits

Self-collect kits available from CHLabs (patient call to collect from CHLabs reception).

Kit contains collection device with swab, transfer pipette & urine collect pottle for self-take vaginal swab or urine sample.

## **Virology Swabs**

For stock orders email labinfo@cdhb.health.nz. For other enquires phone Virology: ext. 80356

- 1. Herpes simplex virus, Varicella zoster Rash or ulcer swabs in viral transport medium.
- 2. Mpox Consult Microbiologist. Take swabs from multiple lesions. Double bag with any other samples taken.
- 3. Syphilis Genital ulcer swab in viral transport medium.
- 4. Adenovirus PCR Eye swab in viral transport medium.
- 5. Enterovirus PCR, Mycoplasma pneumoniae PCR Throat swab in viral transport medium.
- 5. Mumps PCR Buccal, oral, saliva swab in viral transport medium, or CSF neat.
- 6. Respiratory pathogens incl. Flu, RSV, SARS-CoV-2, Bordetella pertussis, and Measles -Nasopharyngeal swab in viral transport medium. Paediatric swabs available on request.







Swab has a breakpoint; break swab carefully into transport medium.



