

This document is valid at the date of issue and will be updated as required within Canterbury Health Laboratories. People outside of Canterbury Health Laboratories should use this document with caution, or contact the laboratory for updated information.

General Chemistry

				Date Revised
Sodium All ages	135	145	mmol/L	Jun-09
Potassium All ages	3.5	5.2	mmol/L	Jun-09
Chloride	95	110	mmol/L	Jun-09
Bicarbonate 0-2yrs >2yr - Adult	19 22	24 32	mmol/L	Dec-21
Anion Gap	10	20	mmol/L	Dec-04
Glucose- Random	3.5	7.7	mmol/L	Jun-09
Glucose - Fasting	3.5	6.0	mmol/L	Jun-09
Urea 1 day - 30 days >1 mth - 2yr >2 yrs - 4yrs >4 yrs - 6 yrs >6 yrs - 10 yrs >10 yrs - 13 yrs >13 - yrs Adult	1.1 0.7 1.1 1.1 1.4 1.8 3.2	6.1 5.0 5.0 5.7 5.7 6.4 7.7	mmol/L	Aug-14
Creatinine 0 - 1 month >1 mth - 2 yr >2 yrs - 4 yrs >4 yrs - 6 yrs >6 yrs - 10 yrs >10 yrs - 15 yrs Male >15yr + Adult Female >15 yr + Adult	20 20 20 25 25 40 50 45	60 50 60 65 70 80 110 90	µmol/L	Aug-14
est GFR	80	120	ml/min/1.73m ²	Mar-05
Osmolality	280	300	mmol/kg	

12/01/2022

Page 1 of 11

General Chemistry

Urate (Uric Acid)			mmol/L	
0-11 yrs	0.08	0.34		Jun-09
Male	0.20	0.42		
Female	0.15	0.36		
Calcium	2.2	2.6	mmol/L	Jun-09
Calcium - Ionised, serum pH7.4	1.13	1.26	mmol/L	Oct-05
Phosphate			mmol/L	Jun-09
0- 30 days	1.3	2.5		
31 days - 1 yr	1.2	2.2		
> 1 yr - 4 yrs	1.1	2.0		
> 4 yrs - 11 yrs	1.0	2.0		
> 11 yrs - 16 yrs	0.9	1.9		
16 yr + Adult	0.8	1.5		
Albumin			g/L	
0-14 days	28	41		Jun-15
>14 d -1 yr	28	45		
>1yr - 8yr	35	45		
>8yr - 15 yr	37	47		
> 15yr - Adult	32	48		Nov-15
Bilirubin			µmol/L	Aug-14
Cord blood	10	40		
Bilirubin				
age same as DOB	0	50		
1 day	0	100		
2 days	0	140		
3 days	0	200		
4-14 days	0	170		
>2 - 6 weeks	3	23		
>6weeks	2	20		
Conjugated Bilirubin	0	5	µmol/L	
Alkaline Phosphatase (ALP)			U/L	
0-1 yr	60	500		Aug-14
>1-11 yrs	80	360		
Male >11-16yrs	90	450		
Female >11-16 yrs	80	400		
Male >16-18 yrs	60	300		
Female >16-18 yrs	50	200		
> 18 + Adult	30	150		

General Chemistry

Gamma glutamyl transferase (GGT)			U/L	Aug-14
0 - 6 mths		<150		
>6mth - 1 yr		<50		
>1 yr - 12 yrs		<30		
>12 yr - Adult Male	10	50		
>12yr- Adult Female	10	35		
Lactate	0.5	2.0	mmol/L	Jul-08
C-Reactive Protein (CRP)		<5	mg/L	Jun-09
β Hydroxybutrate		<0.3	mmol/L	Nov-13
Total Protein			g/L	Aug-14
0 - 2 mth	45	70		
>2mth - 1 yr	50	80		
>1yr - 6 years	55	80		
> 6yrs - 11 yrs	60	83		
> 11 yrs -adult	64	83		
AST (Aspartate aminotransferase)			U/L	Aug-14
0 - 1 month	20	200		
>1 mth - 2 yr	20	80		
>2 years - 5 years	15	60		
> 5 years + adult	10	50		
ALT (Alanine aminotransferase)			U/L	Jun-09
Male	0	40		
Female	0	30		
Magnesium	0.6	1.2	mmol/L	Jun-09
LD (Lactate Dehydrogenase)(L-P)	(M & F)		U/L	Sept-17
0-14d	310	1220		
>14d - 1yr	190	450		
>1yr - 10yr	190	320		
>10yr - 15 yr	160	280		
>15yr - 18 yr	130	250		
>18yr	110	220		

General Chemistry

CK (Creatine Kinase)			U/L	Jul-13
Male	60	220		
Female	30	180		
Amylase (Total)		<100	U/L	Jun-09
Pancreatic Amylase	8	53	U/L	
Lipase	10	70	U/L	Jun-13
Ammonia		<50	µmol/L	May-10
< 30 days				
> 30days, Adult	11	32		
Bile Acids (Non-fasting)		<15	µmol/L	Sept-19

Lipids

Cholesterol	<4.0	mmol/L	Jun-09
Triglyceride	<1.7	mmol/L	Jun-09
HDL Cholesterol	>1.0	mmol/L	Jun-09
LDL Cholesterol (Calculated)	<2.5	mmol/L	Jun-09
Cholesterol/HDL Cholesterol ratio	<4.5	ratio	Jun-09

Optimal lipid levels are suggested by the New Zealand Guidelines Group <http://www.health.govt.nz/publication/new-zealand-primary-care-handbook-2012> though an individualised multifactorial approach is preferred.

Iron Studies

Iron			µmol/L	Aug-14
0 - 2 months	10	31		
>2 mths - 1 year	4	27		
>1 - 3 yrs	5	23		
>3 - 10 years	6	25		
>10 - 18 years	8	32		
>18yr + adult	10	30		
Transferrin			g/L	Aug-14
0 - 2 months	1.3	2.3		
>2 mths - 3 years	1.6	3.3		
>3yr + Adult	2.0	3.5		
Transferrin Saturation (Iron x 3.98/ Transferrin)	16	45	%	
Ferritin			µg/L	Sep-21
0 - 15 years	15	150		
Male >15 - 30 years	20	350		
Male >30 years	20	400		
Female >15 - 50 years	20	150		
Female >50 years	20	300		
Vitamin B12	105	675	pmol/L	Sept-21
Folate - Serum		>7	nmol/L	Sept-21

Cardiac Markers

BNP		<30	pmol/L	Dec-05
CKMB		<5.0	µg/L	Oct-05
Myoglobin				Oct-05
Male		<155	µg/L	
Female		<106	µg/L	
hs TNI				
Male	0	20	ng/L	Sept-21
Female	0	10	ng/L	
High Sensitive C-Reactive Protein (hsCRP)		<3	mg/L	Nov-17
<hr/>				
Foetal Fibronectin		<50 Negative	ng/mL	May-17
		≥50 Positive		
Pregnancy β-HCG		<5	Negative	U/L
		25	Equivocal	
		>25	Positive	

Tumour Markers

Total β-HCG		<5	U/L	Aug-18
AFP	0	16	ug/L	Aug-18
CEA	0	5	ug/L	Aug-18
CA 125	0	35	kIU/L	Aug-18
CA 15-3	0	24	kIU/L	Sept-21
CA 19-9	0	35	kIU/L	Sept-21
PSA	0	4.0	ug/L	

Thyroid Function Tests

T3 (Free)				
0 - 1yr	4.3	6.9	pmol/L	Sept-21
>1 -15 years	4.0	6.2		
>15 - Adult	3.3	6.8		
T4 (Free)				
0- 20 days	13	58	pmol/L	Dec-21
>20dy - 3 years	9	18		
>3yrs - Adult	8	20		
TSH				
1 - 4 days	0.8	20.0	mU/L	Sept-21
>4 dy - 1 mth	0.8	12.0		
>1 mth - 6 mths	0.8	5.9.0		
> 6 mths + Adult	0.4	5.3.0		
Anti-TPO	<10		U/mL	Sept-21

(Anti – TG removed Feb -20)

Vitamin D

Optimum Target Range for bone health	50	150	nmol/L	Feb-19
--	----	-----	--------	--------

25-OH-Vitamin D levels < 25 nmol/L are indicative of moderate to severe vitamin D deficiency.

25-OH-Vitamin D levels 20 – 50 nmol/L are indicative of mild vitamin D deficiency

25-OH-Vitamin D levels > 200 nmol/L are found in vitamin D toxicity.

Serum Indices(c16000)

Haemolysis	0	34	Specimen not haemolysed
	35	99	Specimen slightly haemolysed
	100	249	Specimen haemolysed
		≥250	Specimen grossly haemolysed
Icterus Index		>40	Specimen icteric
Lipaemic Index		>3.3	Specimen lipaemic

Drug Monitoring

Digoxin	0.6	2.0	nmol/L	Jun-09
Salicylate	1.5	1.8	mmol/L	
Paracetamol	umol/L	as per latest edition of Hospital Health Pathways		
Ethanol	<3	Undetected	mmol/L	
Gentamicin (trough)	<0.5		mg/L	Aug-09
AUC	70	100	mg/L/hr	
Tobramycin (trough)	<0.5		mg/L	Aug-09
AUC	70	100	mg/L/hr	
Vancomycin	5	10	mg/L	Aug-09
Carbamezapine	16	50	µmol/l	Aug-09
Unbound	4	10	µmol/L	
Lamotrigine	12	55	µmol/L	Dec-15
Phenobarbitone	65	130	µmol/l	Aug-09
Phenytoin	40	80	µmol/L	Aug-09
Unbound	4	9	µmol/L	
Theophylline	55	110	µmol/L	Aug-09
Valproate Adult	350	700	umol/L	Aug-09

Blood Gases – Whole Blood

	Arterial	Venous	Units	
pH	7.35 - 7.45	7.30 - 7.40		Jan-20
pCO ₂	35 - 45	40 - 50	mmHg	Jan-20
pO ₂	75 - 105		mmHg	Jan-20
p50		22 - 27	mmHg	
HCO ₃	22 - 28	22 - 32	mmol/L	Jan-20
Base Excess	±3	±3	mmol/L	
Total CO ₂	24.0 - 30.0	25.0 - 33.0	mmol/L	
O ₂ Saturation	>97		%	
Lactate	< 2.0	< 2.0	mmol/L	Jan-20
Ionised Calcium	1.15 - 1.30	1.15 - 1.30	mmol/L	Jan-20
Ionised Calcium pH corrected 7.4	1.13 - 1.26	1.13 - 1.26	mmol/L	May-10
tHb	Refer Core Haematology			May-10
Oxy Hb	94 - 97		%	
Carboxy Hb non-smokers	0 - 1.5	0 - 1.5	%	
smokers	<9.0	<9.0	%	
MetHb	0 - 1.5	0 - 1.5	%	
DeOxy Hb	0 - 5.0		%	

Blood Gases – Cord Blood

	Cord Arterial	Cord Venous	Units	
pH	7.10 - 7.40	7.22 - 7.44		Jan-20
pCO ₂	42 - 74	28 - 53	mmHg	Jan-20
pO ₂	6 - 28	16 - 40	mmHg	Jan-20
HCO ₃	18 - 28	17 - 26	mmol/L	Jan-20
Base Excess	-12 to +2	-12 to +2	mmol/L	Jan-20
Lactate	< 6.1	<6.1	mmol/L	Jan-20

CSF Chemistry

Appearance	Clear & Colourless		
Protein	0.15	0.4	g/L
Glucose	2.8	4.4	mmol/L
Creatinine	20	50	µmol/L
LD			U/L
Lactate	1.1	2.3	mmol/L
Xanthochromia	Not detected		

Urine Chemistry

Sodium	100	250	mmol/24hr	
Potassium	35	100	mmol/24hr	
Chloride	110	250	mmol/24hr	
Calcium	2.5	7.5	mmol/24hr	
Calcium/Creatinine Ratio	0.06	0.45	mole/ratio	
Phosphate	16	48	mmol/24hr	
TRP (Tubular Reabsorption of Phosphate)	82	97	%	2003
Glucose		<1	mmol/24hr	
Creatinine			mmol/24hr	Aug-14
0 - 6 months	0.2	1.6		
>6mths - 2 years	0.2	2.0		
>2 - 6 years	0.5	6.0		
>6 - 10 years	1.5	12.6		
>10 - 15 years	2.0	20.0		
Male >15 years	9.0	18.0		
Female >15years	7.0	12.5		
Osmolality				
0 - 6 years	50	1200	mm/Kg	
>6 yrs + Adult	300	1200		
Protein	0	0.14	g/24Hr	
UPRO/CRN ratio		<22.7	g/mol	
Albumin		<30	mg/L	
UALB/CRN ratio				
Male		<2.5	g/mol	
Female		<3.5	g/mol	
Urea	165	580	mmol/24hr	
Amylase		<35	U/mmol CRN	
Creatinine Clearance		>1.5	ml/sec	

Copies of this Document are kept as follows:

- 1 APS Operations Manual
- 2 Manual Test Procedures
- 3 Checkout Bench
- 4 CE Bench
- 5 Section manual
Clinical
- 6 Haematology
Mark Warner -
- 7 CCST
Research Team
- 8 CHL

E-mail copies to:

Biochemistry Peer Review Group

**Reference Range
updates required:**

- Instrument
- 1 Manager
 - 2 POCT devices, ABL90, iStat